



Lichtmangel ausgleichen: Wohlfühlen im Winter

Im Herbst und Winter kann es zu Lichtmangelercheinungen kommen. Oft werden diese Störungen des Hormonhaushalts als Winterdepression bezeichnet. Symptome sind:

- Unausgeglichenheit
- gedrückte Stimmung
- Energie- und Antriebslosigkeit
- allgemeines Unwohlsein

Künstliches Licht kann dies nicht ausgleichen. Das spezielle, helle Licht der Tageslichtlampe kann als Lichtdusche - ohne Medikamente - auf den Hormonhaushalt einwirken und vorbeugend oder behandelnd eingesetzt werden.

Balancing shortage of light: Well-Being in winter

In winter, there may be light deficiency (disturbance of the hormone-budget):

- imbalance
- depressed mood
- less energy and motivation
- general malaise

Artificial light can not compensate.

The special, bright light of daylight lamp can influence - without medicine - the hormone budget and preventive or treated used.

german|engineering

TL 40 Tageslichtlampe

Lichtdusche mit Tageslicht
- zur Anwendung bei Lichtmangelercheinungen

Ideal für den Schreibtisch geeignet

Simulation von Sonnenlicht:

Lichtstärke ca. 10.000 Lux
(bei 15 - 20 cm Abstand)

Besonders helle und gleichmäßige Ausleuchtung

Beleuchtungsfläche ca. 21 x 20 cm

Flackerfreie Zündung - sofort angenehmes Licht
ohne Stroboskopeffekte

Flimmerfrei und UV-frei

Energiesparendes Leuchtmittel

Mit Aufsteller

Zum Aufhängen geeignet

Komfortable Ein-Knopf-Bedienung

Zertifiziertes Medizinprodukt

1 x 45 Watt Röhre

Wirkungsgrad des Vorschaltgerätes: ~ 96%

Wirkungsgrad der Lampe: 67-71 lm/W

Farbtemperatur der Leuchtröhren: 6.500 Kelvin

Produktmaße: 20 x 7 x 27 cm

Produktgewicht: ca. 800 g

3 Jahre Garantie

VE: 4 / Transportkarton: -

EAN-Nr.: 4211125 60820 0

Art.-Nr.: 608.20

TL 40 Brightlight

Light box with daylight
- For use in case of light deprivation

Ideal suitable for your writing desk

Simulation of sunlight:

Intensity of light approx. 10.000 Lux
(at a distance of 15 - 20 cm)

Exceptionally bright and even illumination

Illumination surface approx. 21 x 20 cm

Flicker-free starting - immediately pleasant light
without stroboscopic effects

Flicker-free and UV-free

Energy-saving light source

With flap

Suitable for hanging

Convenient one-button operation

Certified medical product

1 x 45 Watt tube

Degree of efficiency of the ballast: ~ 96%

Source efficiency of the bulb: 67-71 lm/W

Color temperature of tubes: 6.500 Kelvin

Product measurements: 20 x 7 x 27 cm

Product weight: approx. 800 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 60820 0

Item no.: 608.20